

Simply Cloning

A step-by-step video guide for making DNA constructs



Cloning Strategies

www.cloningstrategies.com

Molecular Cloning Protocols

1 | Planning

Insert

2 | PCR

4 | Gel Purification

5 | Restriction Digest

Vector

3 | Restriction Digest

4 | Gel Purification

6 | Ligation

7 | Transformation

8 | Plasmid Miniprep

9 | Verification: Restriction, Sequencing

10 | Depositing Clones in a Library

1 | Planning

- o Create experiment outline in PowerPoint. Include:
 - Name of starting vector
 - Insert template
 - Restriction enzymes
 - Primers
- o Build vector maps
- o Design primers:
 - Select first/last 20 nucleotides of the insert sequence
 - Add a restriction site at the 5' end of the forward primer and at the 3' end of the reverse primer
 - Add six random nucleotides to enable nuclease cleavage close to the end of DNA fragment
 - For the forward primer: insert one or two nucleotides between the restriction site and the insert sequence to maintain the reading frame of your fusion protein
 - For the reverse primer: add a stop codon if necessary
 - For the reverse primer: transform to complimentary sequence

Forward Random *Bam*HI Frame Insert sequence
primer: 5' CGATGA GGATCC GT GGTAGCAAGGGCGAGGAGC 3'

Reverse Insert sequence Stop *Spe*I Random
primer: 3' GTACCTGCTCGACATGTTT TGA ACTAGT ACATAG 5'

Reverse
complim.: 5' CTATGTACTAGTTCAGAACATGTCGAGCAGGTAC 3'

2 | PCR

- o PCR mix:
 - 37.5 μ l ddH₂O
 - 5 μ l 10X PCR buffer
 - 4 μ l dNTPs (10mM each)
 - 1 μ l Primer (forward), 10 μ M
 - 1 μ l Primer (reverse), 10 μ M
 - 1 μ l DNA template
 - 0.5 μ l Polymerase (1.25 units)
- Note: If you are using a plasmid miniprep as your DNA template, dilute it 50 times and take 1 μ l of diluted solution

- o PCR program:
 1. 94 °C – 2 min
 2. 94 °C – 30 sec
 3. 55 °C – 30 sec
 4. 72 °C – 1 min (or 1 min/kb)
 5. Go back to 2. Repeat 29 times
 6. 72 °C – 5 min
 7. 4 °C – forever (insert 0 on most thermal cyclers)

3 | Vector Restriction Digest

- o Restriction mix:
 - 39 μ l ddH₂O
 - 5 μ l 10X enzyme buffer
 - 5 μ l DNA (approximately 1 μ g)
 - 0.5 μ l restriction enzyme A (5 units)
 - 0.5 μ l restriction enzyme B (5 units)
- o Mix up by pipetting up and down
- o Incubate at 37 °C for 45 min

4 | Gel Purification

- o Prepare 1% agarose gel:
 - Add 1 g of agarose to 100 ml of TEA buffer
 - Bring to boil in a microwave, make sure the agarose completely dissolved
 - Cool at room temperature for 10 min
 - Add 2 µl ethidium bromide, pour into a tray
 - Leave for 30 min to solidify
- o Add 5 µl of loading buffer (60% glycerol in TAE, trace amount of bromophenol blue) to 50 µl of your restriction digest/PCR
- o Load gel:
 - Molecular weight marker
 - 1 µl of undigested vector (mixed with loading dye)
 - Vector restriction digest
 - Insert restriction digest
- o Run gel until bromophenol blue dye reaches the end
- o Take gel picture
- o Cut out desired DNA fragments
- o Purify with a gel purification kit. Elute with 50 µl of ddH₂O

7 | Transformation and Plating

- o Thaw chemically competent cells on ice for 15 min
- o Add ligation mix to the competent cells, mix by tapping gently
- o Leave the cells on ice for 30 minutes
- o Heat-shock cells at 42 °C for 60 seconds
- o Put cells back on ice for 5 minutes
- o Add 4 volumes of LB medium without antibiotics (for example, to 200 µl of competent cells add 800 µl of LB)
- o Shake for 45 min at 37 °C for transformed cells to develop antibiotic resistance
- o Spin down at 13000 rpm for 1 min, remove 800 µl of supernatant, resuspend pelleted cells in the remaining LB
- o Plate on LB agar plate with antibiotics
- o Incubate the plate at 37 °C overnight

9 | Clone Verification by Restriction Digest

- o Restriction mix:
 - 16.6 µl H₂O
 - 1 µl miniprep DNA
 - 2 µl 10X restriction buffer
 - 0.2 µl restriction enzyme A (2 units)
 - 0.2 µl restriction enzyme B (2 units)
- o Incubate at 37 °C for 45 min
- o Run on an agarose gel, take a picture

5 | Insert Restriction Digest

- o Restriction mix:
 - 44 µl of gel purified insert from step I-2
 - 5 µl 10X Buffer
 - 0.5 µl restriction enzyme A (5 units)
 - 0.5 µl restriction enzyme B (5 units)
- o Mix up by pipetting up and down
- o Incubate at 37 °C for 45 min
- o Deactivate the enzymes:
 - A. Clean on a GFX column (enzyme removal kit)
 - B. Check the enzymes you are using if they can be heat deactivated. For most enzymes, incubate 30 min at 65 °C

6 | Ligation

- o Ligation mix:
 - 4 µl H₂O
 - 1 µl 10X ligation buffer
 - 1 µl of plasmid digested with restriction enzymes A and B
 - 3 µl of insert digested with restriction enzymes A and B
 - For the ligation control, use 3 µl of water instead of insert
 - 1 µl DNA ligase
- o Incubate at 16 °C for 15 min - 24 hours
Note: longer incubation times increase efficiency

8 | Plasmid Miniprep

- o Prepare numbered cell culture tubes with 4 ml of liquid LB medium and appropriate antibiotics
- o Pick colonies from the transformation plate with autoclaved toothpicks or pipette tips and inoculate the corresponding media in the numbered tubes
- o Incubate the tubes in a shaker at 37 °C overnight
- o On the next day, save 1 ml of culture in a separate eppendorf tube (for clone library)
- o Spin down the remaining 3 ml of culture and proceed with your plasmid miniprep kit

10 | Depositing Clones in a Library

- o Create an Excel or Word file with library information
Include:
 - Clone number and name
 - Clone description
 - Primers used
 - Date created
 - Anything else you find relevant
- o For each clone, prepare two tubes with clone name and number
- o Add 0.9 ml of overnight E. coli culture and 0.1 ml of DMSO to one of the tubes, store at -80 °C
- o Put plasmid DNA into another tube, store at -20 °C